

F.A.N. PAD-GL™

Formaldehyde Absorbent Neutralizer Pad

“ULTRA-SHIPPER”

The Revolutionary Transport Insert Pad

**Ship & Receive Specimen Containers
With Confidence**

**No More Harmful Vapors
From Leaking or Damaged Containers
While Transporting Formalin Pre-filled Specimen Containers.**



**Insert The New F.A.N. Pad GL™ “ULTRA” SHIPPER Pad
and Specimen Container into the Bio-Hazard Transport Bag**
Each F.A.N. Pad-GL™ “ULTRA” 4” x 4” Shipper Pad will absorb and neutralize 20-ml of Formalin

If there are No leaks during transport,
Place a tissue cassette on the pad and
Pour off specimen into the cassette.
Pad will absorb & neutralize formalin.

OR

Place Larger Specimen on the Pad
for Dissecting the Tissue
Excellent Knife Suppression,
Increases life of Dissection Blades.

▶ Reduces Formalin Vapor Exposure

▶ Reduces Transport Liability

▶ Non-Hazardous Pad Disposal.

▶ Reduces Risk of Lost Specimen

▶ **No More “NASTY POOF”, Opening Specimen Bags.**

Product # FAN- 404 200 - 4” x 4” Pads/ Carton

American Bio-Safety, Inc.

4322 Anthony Court, # 5 Rocklin, CA. 95677

Bus: 916.652.8021 Email: absafety@sbcglobal.net Web Site: www.AmericanBioSafety.com Fax: 916.652.8020